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## Anexo

**BIOLABO REAGENTS**  
www.biolabo.fr

**MANUFACTURER:**  
**BIOLABO SA,**  
02160, Maizy, France

# CALCIUM

## Arsenazo III method

Reagent for quantitative determination of calcium in human serum and plasma, or urines.

REF 90004 R1 2 x 125 mL R2 1 x 10 mL

### TECHNICAL SUPPORT AND ORDERS

Tel : (33) 03 23 25 15 50

Fax : (33) 03 23 256 256

IVD IN VITRO DIAGNOSTIC USE

### CLINICAL SIGNIFICANCE (1) (2)

Calcium fulfills a variety of roles in human physiology, not only as a structural factor in bones and teeth, but also in normal neuromuscular function and clotting of blood.

The level of serum calcium may be affected by intestinal malabsorption, by alterations in plasma proteins level, especially albumin, which should be measured concurrently with calcium.

Hypercalcemia is found in hyperparathyroidism, multiple myeloma, bone and parathyroidal neoplasms and in states with bones demineralisation.

Hypocalcemia is encountered in hypoparathyroidism and in several cases of necrosis and acute pancreatitis.

### PRINCIPLE (4)

At mildly acidic pH, metallo-chromogen Arsenazo III combines with calcium to form a coloured complex which absorbance measured at 650 nm (640-660) is proportional to the amount of calcium in the specimen.

### REAGENTS

Vial R1

#### ARSENAZO III REAGENT

Imidazol buffer pH 6.8 at 25°C	> 90	mmol/L
Arsenazo III	> 0.18	mmol/L
Surfactant	0.1	%
Preservative		

Vial R2

#### STANDARD

Calcium 10 mg/dL (2.5 mmol/L)

### SAFETY CAUTIONS

BIOLABO reagents are designated for professional, in vitro diagnostic use.

- Use adequate protections (overall, gloves, glasses).
- Do not pipette by mouth.
- In case of contact with skin and eyes, thoroughly wash affected areas with plenty of water and seek medical advice.
- Reagents contain sodium azide (concentration < 0.1%) which may react with copper and lead plumbing. Flush with plenty of water when disposing.
- Material Safety Data Sheet is available upon request.
- Waste disposal : Respect legislation in force in the country.

All specimens should be handled as potentially infectious, in accordance with good laboratory practices using appropriate precautions. Respect legislation in force in the country.

### REAGENT PREPARATION

Reagents are ready for use.

### STABILITY AND STORAGE

Store at 18-25°C, away from light.

- Unopened :  
Reagents are stable until expiry date stated on the label.
- Once opened :  
Reagent R1 is stable at least for 3 months when free from contamination.

Standard stability (vial R2) : Several weeks once opened (transfer requested quantity, recap and store at 18-25°C).

Discard any reagent if cloudy or if reagent blank at 650 nm > 0.400.  
This kit can travel at room temperature.

### SPECIMEN COLLECTION AND HANDLING (1) (2)

Serum or heparinised plasma :

Do not use citrate, oxalate or EDTA. Blood obtained on fasting patient with minimal venous occlusion and without exercise or after restoring circulation at least for 1 minute.

24 h Urines :

Acidify after collection with 20 to 30 mL HCl 6 N to dissolve calcium salts.

Dilute (1 + 2) with distilled water before performing the test.

Total calcium is stable in serum for :

- at least 7 days at 2-8°C.
- 6 months at -20°C.

Long-term freezing may lead to associated evaporation, lyophilisation or coprecipitation with fibrin (i.e. heparinised plasma) or lipids.

### INTERFERENCES (3)

Hemolysis, icterus, lipemia, paraproteins and magnesium : Perform bichromatic or multi-wavelength analysis or specimen blank to reduce positive or negative interferences.

Handle with care specimens, calibrators and controls to avoid contamination by environmental calcium. Use disposable tubes and cuvettes and clean glassware with HCl 1N, well rinse with demineralised water.

Plastic and glass containers may adsorb calcium during storage, especially with diluted solution.

For a more comprehensive review of factors affecting this assay refer to the publication of Young D.S.

Version : AT 90004 21 02 2005

## MATERIAL REQUIRED BUT NOT PROVIDED

1. Basic medical analysis laboratory equipment.
2. Normal and pathological control sera.

## CALIBRATION

- Kit Standard (vial R2) or BIOLABO-Multicalibrator, REF 95015.
- Or any calibrator traceable to a reference method or material.

The calibration frequency depends on proper instrument functions and on the preservation of the reagent.

It is recommended to calibrate in the following cases :

1. When changing batch of reagent.
2. After maintenance operations on the instrument .
3. If control values are out of range, even after using a new vial of fresh serum.

## QUALITY CONTROL

- BIOLABO EXATROL-N (normal values), REF 95010.
- BIOLABO EXATROL-P (pathological values), REF 95011.
- Assayed control sera referring to the same method.
- External quality control program.

It is recommended to control in the following cases :

- At least once a run.
- At least once within 24 hours.
- When changing vial of reagent.
- After maintenance operations on the instrument.

If control is out of range, apply following actions :

1. Repeat the test with the same control.
2. If control is still out of range, prepare a fresh control serum and repeat the test.
3. If control is still out of range, use a new vial of calibrator or a fresh calibrator and repeat the test.
4. If control is still out of range, calibrate with a new vial of reagent.
5. If control is still out of range, please contact BIOLABO technical support or your local Agent.

## EXPECTED VALUES (2)

### TOTAL CALCIUM in serum :

Population	mg/dL	mmol/L
Premature	6.2-11.0	[1.55-2.75]
0-10 days	7.6-10.4	[1.90-2.60]
10 days –24 months	9.0-11.0	[2.25-2.75]
24 months –12 years	8.8-10.8	[2.20-2.70]
12 years -18 years	8.4-10.2	[2.10-2.55]
18-60 years	8.6-10.0	[2.15-2.50]
60-90 years	8.8-10.2	[2.20-2.55]
> 90 years	8.2-9.6	[2.05-2.40]

**TOTAL CALCIUM in 24 h Urines :** < 300 mg/24 h (< 7.5 mmol/24 h).

Each laboratory should establish its own normal ranges for the population that it serves.

## PERFORMANCES CHARACTERISTICS

Within run n = 20	Normal Level	High Level	Between run n = 20	Normal Level :	High Level
Mean mg/dL	8.82	14.02	Mean mg/dL	8.96	14.14
S.D. mg/dL	0.04	0.06	S.D. mg/dL	0.11	0.20
C.V. %	0.5	0.4	C.V. %	1.26	1.42

Detection limit : approximately 0.21 mg/dL

Sensitivity for 1 mg/dL : 0.070 Abs. at 650 nm.

Comparison with commercially available reagent :

$$y = 1.0259 x - 0.1811 \quad r = 0.9943$$

## LINEARITY

The assay is linear up to 15.0 mg/dL (3.75 mmol/L). Above, dilute the specimen with demineralised water and re-assay taking into account the dilution factor. Linearity limit depends on specimen/reagent ratio.

## MANUAL PROCEDURE

Temperature should be held constant as the absorbance of the dye is temperature sensitive.

Pipette into well identified test tubes :	Blank	Standard	Assay
Reagent	1 mL	1 mL	1 mL
Distilled water	20 µL		
Standard		20 µL	
Specimen			20 µL

Mix well. Let stand for 1 minute at room temperature.  
Read absorbances at 650 nm (640-660) against reagent blank.  
The coloration is stable for 1 hour away from light

### Notes :

- ✓ Bichromatic analysis : the 2<sup>nd</sup> wavelength is 700 nm.
- ✓ Specific procedures are available upon request for automated instruments. Please contact BIOLABO technical support.

## CALCULATION

Calculate the result as follows :

$$\text{Serum or plasma :} \quad \text{Result} = \frac{\text{Abs(Assay)}}{\text{Abs(Standard)}} \times \text{Standard concentration}$$

Urines diluted (1 + 2) : multiply the above result by dilution factor 3.

## REFERENCES

- (1) TIETZ N.W. *Text book of clinical chemistry*, 3<sup>rd</sup> Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p. 1395-1406, p.1435-1439.
- (2) *Clinical Guide to Laboratory Test*, 3<sup>rd</sup> Ed., N.W. TIETZ (1995) p. 102-104.
- (3) YOUNG D.S., *Effect of Drugs on Clinical laboratory Tests*, 4<sup>th</sup> Ed. (1995) p. 3-115 à 3-125
- (4) BAUER J. P., *Affinity and stoichiometry of calcium binding Arsenazo III*, *Anal. Biol. Chem.*(1981), 110, p.61-72

Made in France

Version : AT 90004 21 02 2005

Manufacturer

Use by

IVD In vitro diagnostic

Temperature limitation

REF Catalogue number

See insert

LOT

Batch number